

### ***Rule 1***

**Follow the directions in your course reader, of your teaching assistant and of your instructor. They are usually much more experienced doing chemistry.**

### ***Rule 2***

**Be prepared when entering the lab. It is for your own and everybody else's safety. This will also increase your chances to complete the experiments in a timely fashion.**

### ***Rule 3***

**When in doubt, ask. This will make it more likely that you isolate the product and keep your health.**

### ***Rule 4***

**Keep all reagent and solvent bottles closed when not in use! This will keep the reagents in better condition for everybody.**

### ***Rule 5***

**Always inspect your glassware before using it. Initial cracks can lead to breakage during heating or implosions upon evacuation.**

### ***Rule 6***

**Lubricate ground-glass joints lightly when a tight seal is needed (exclusion of moisture or vacuum).**

### ***Rule 7***

**When running a reaction, a stir bar or spin vane (point down) has to be added which has to be spinning, of course.**

### ***Rule 8***

**The O-ring is placed between the compression cap and the conical vial.**



### ***Rule 9***

**Never heat anything up in a closed system. It will explode eventually because the vapor pressure of every compound increases upon heating, the question is only how quickly.**

### ***Rule 10***

**Do not fill the vessel more than half when refluxing or distilling. There has to be room for bubbles and mixing.**

### ***Rule 11***

**Reflux means to boil a liquid in a container equipped with a cooled condenser so that the vapors continuously condense for re-boiling.**

### ***Rule 12***

**An ice-bath consists mainly of water, which has some ice added to reduce temperature.**

### ***Rule 13***

**The container and the heat (or cooling) source have to have good contact for heat transfer to occur.**

### ***Rule 14***

**Acid-base reactions are faster than nucleophilic additions. Thus, protic solvent cannot be used with strong nucleophiles, which are often also strong bases.**

### ***Rule 15***

**Hydrogen bonds are stronger than dipole-dipole interactions and van der Waals interactions but they are still much weaker than covalent bonds.**

### ***Rule 16***

**The stronger the intermolecular bonds are the higher the melting point and the boiling point of a compound.**

***Rule 17***

Symmetric molecules have higher melting points than asymmetric molecules even if they have smaller dipole moments.

***Rule 18***

Intramolecular reactions will be favored if five-membered or six-membered rings are formed.

***Rule 19***

An increase of ten degrees Celsius causes a doubling of the rate of the reaction.

***Rule 20***

Never dispose of any layer, until you are sure that you will never need it again.

***Rule 21***

Use a short-stem funnel when pouring the solution into the separatory funnel, which should not be filled more than two-thirds.

***Rule 22***

The compound should exhibit a high solubility in the solvent used for extraction for the extraction to work well.

***Rule 23***

About 10-20% of the volume of the solution that is being extracted or washed should be used for each extraction.

***Rule 24***

Multiple extractions (2-3) with small amounts are more effective than one extraction with the same total volume. The first extract contains most of the product.

***Rule 25***

Acidic compounds are extracted with bases, basic compounds with acids. Water is extracted with saturated sodium chloride solution.

***Rule 26***

Vent the separatory funnel frequently, particularly when using sodium bicarbonate or low boiling solvents for the extraction.

***Rule 27***

When separating layers during extraction, the bottom layer is always removed first independent if this is the layer of interest or not.

***Rule 28***

Drying agents have to be used sparingly because the more polar a compound is the higher the affinity towards the drying agent will be.

***Rule 29***

The drying agent has to be removed from a solution before proceeding to the next step because the drying process is reversible at elevated temperatures.

***Rule 30***

In a distillation or a reflux setup the vent has to be placed after (or on top of) the condenser to ensure that the condenser works properly.

***Rule 31***

Excessive clamping usually leads to breakage because of strain. In most cases, two fixed points in the proper location will be sufficient to support the setup.

***Rule 32***

Never distill to dryness because this will not lead to any separation. It can cause an explosion because many peroxides have higher boiling points than the actual compound.

***Rule 33***

Compounds exhibiting a boiling point above 150 °C have to be distilled *in vacuo* to prevent decomposition.

***Rule 34***

Boiling stones cannot be used in recrystallization and during vacuum distillations.

***Rule 35***

When recrystallizing the target compound and the solvent have to have different polarities, while the impurity and the solvent should have the similar polarity.

***Rule 36***

The solution has to be brought to a boil during the process to maximize the solubility of the compound and minimize the amount of the solvent required to dissolve the sample.

***Rule 37***

The entire crude has to be dissolved during the process. Undissolved materials have to be removed by filtration or decanting.

***Rule 38***

A slow crystallization yields the best product in terms of purity and crystallinity because this will allow the system time to be more selective.

***Rule 39***

Supersaturation can be overcome by scratching on the inside of the vessel or by seeding.

***Rule 40***

When using a mixed solvent system, the solvents have to have different polarities but similar boiling points.

***Rule 41***

Silica and alumina are polar stationary phases. Thus, polar compounds adsorb stronger than non-polar compound on these stationary phases.

***Rule 42***

On polar stationary phases, polar solvents will cause compounds to migrate more than non-polar solvents if the compound dissolves in the solvent.

***Rule 43***

A capillary spotter has to be used to apply compounds to the TLC plate.

***Rule 44***

Samples for gas chromatography have to be much diluted because a capillary column is used.

***Rule 45***

When visualizing spots on the TLC plate, do not look into the UV-light. It will damage your eyes!

***Rule 46***

Infrared spectra that are acquired using an ATR setup have to be corrected because the penetration depth of the infrared beam is highly wavelength dependent.

***Rule 47***

Polar groups give rise to strong peaks in the infrared spectrum because the dipole moment changes more during the motion.

***Rule 48***

Samples that are submitted for NMR analysis have to be placed in a tube of appropriate length and have to be dissolved in one deuterated solvent only.

***Rule 49***

**The entire sample has to dissolve in the solvent in order to be visible in the NMR or UV-Vis spectrum.**

***Rule 50***

**The solvent and the cuvette have to be transparent in the range to be measured.**

***Rule 51***

***Rule 52***

***Rule 53***

***Rule 54***

***Rule 55***

***Rule 56***